

## Anti-CDK9 Rabbit Monoclonal Antibody

Catalog Number: M00794


### Overview

Product Name	Anti-CDK9 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-CDK9 Rabbit Monoclonal Antibody catalog # M00794. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal GIH-3
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50750


### Technical Details

Immunogen	A synthesized peptide derived from human CDK9
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used:</p> <p>WB 1:1000-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:40</p>


## Anti-CDK9 Rabbit Monoclonal Antibody (M00794) Images




Immunofluorescent analysis using the Antibody at 1:50 dilution.



Western blot analysis of CDK9 expression in HeLa cell lysate.




Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunohistochemical analysis of paraffin-embedded human cervix cancer, using CDK9 Antibody.

Immunofluorescent analysis using the Antibody at 1:500 dilution.



## Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDK9 Rabbit Monoclonal Antibody