

Anti-Csk monoclonal Antibody

Catalog Number: M00799-1


Overview

Product Name	Anti-Csk monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Csk monoclonal Antibody catalog # M00799-1. Tested in WB,ICC applications. This antibody reacts with Human,Rat.
Application	ICC, WB
Clonality	Monoclonal 5C10
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P41240


Technical Details

Immunogen	Recombinant full length Human Csk.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 ICC: 1:50

Anti-Csk monoclonal Antibody (M00799-1) Images



Western blot detection of CSK antibody in C6 and Jurkat cell lysates using CSK antibody (1:1000 diluted) .



Immunocytochemistry staining of C6 cells fixed by anhydrous methanol for 2 h at -20°C and using anti-CSK antibody (dilution 1:50).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Csk monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.