

Anti-IL-4Ra / CD124 Reference Antibody (Manfidokimab)

Catalog Number: M00807-2


Overview

Product Name	Anti-IL-4Ra / CD124 Reference Antibody (Manfidokimab)
Reactive Species	Human
Description	Boster Bio Anti-IL-4Ra / CD124 Reference Antibody (Manfidokimab) (Catalog # M00807-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.457EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P24394


Technical Details

Isotype	IgG4SP
Form	Liquid

Anti-IL-4Ra / CD124 Reference Antibody (Manfidokimab) (M00807-2) Images



Anti-IL-4Ra/CD124 Reference Antibody (Manfidokimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-4Ra/CD124 Reference Antibody (Manfidokimab) is more than 99.14%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-4Ra / CD124 Reference Antibody (Manfidokimab)

For Research Use Only. Not for use in diagnostic procedures.