

Anti-PLZF ZBTB16 Monoclonal Antibody

Catalog Number: M00817

Overview

Product Name	Anti-PLZF ZBTB16 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-PLZF ZBTB16 Monoclonal Antibody catalog # M00817. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 5B3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q05516

Technical Details

Immunogen	Purified recombinant fragment of human PLZF expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-PLZF ZBTB16 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.