

## Anti-PLZF ZBTB16 Monoclonal Antibody

Catalog Number: M00817


### Overview

Product Name	Anti-PLZF ZBTB16 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-PLZF ZBTB16 Monoclonal Antibody catalog # M00817. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 5B3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q05516


### Technical Details

Immunogen	Purified recombinant fragment of human PLZF expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000</p> <p>For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a></p>

## Anti-PLZF ZBTB16 Monoclonal Antibody (M00817) Images



Immunofluorescence analysis of HeLa cells using PLZF monoclonal antibody (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Western blot analysis using PLZF monoclonal antibody against HeLa (1) cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLZF ZBTB16 Monoclonal Antibody