

## Anti-PLA2G4A (Phospho-S505) Antibody (Monoclonal, 32P85)

Catalog Number: M00854

### About PLA2G4A

PLA2G4A (Phospholipase A2, Group IVA), is an enzyme that in humans is encoded by the PLA2G4A gene. Tay et al. (1995) mapped the PLA2G4A gene to rat chromosome 13 by PCR-based intercross genotyping and to human 1q25 by fluorescence in situ hybridization. By site-directed mutagenesis and biochemical analysis of the recombinant protein, Sharp et al. (1994) determined that ser228 participates in the catalytic mechanism of cPLA2 and that both the phospholipase A2 and the lysophospholipase activities are catalyzed by the same active site residue(s). PLA2G4A, the cytosolic phospholipase A2, appears to subserve transmembrane signaling responses to extracellular ligands (Skorecki, 1995).

### Overview

Product Name	Anti-PLA2G4A (Phospho-S505) Antibody (Monoclonal, 32P85)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-PLA2G4A (Phospho-S505) Antibody (Monoclonal, 32P85) catalog # M00854. Tested in WB, DotBlot applications. This antibody reacts with Human, Mouse.
Application	Dot blot, WB
Clonality	Monoclonal 32P85
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P47712

### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser505 of human cPLA2 aa 493-538.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLA2G4A (Phospho-S505) Antibody (Monoclonal, 32P85)

For Research Use Only. Not for use in diagnostic procedures.