

Anti-SPARC Mouse Monoclonal Antibody [Clone ID: OTI2B1]

Catalog Number: M00862


Overview

Product Name	Anti-SPARC Mouse Monoclonal Antibody [Clone ID: OTI2B1]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio SPARC mouse monoclonal antibody, clone OTI2B1 (formerly 2B1). Catalog# M00862. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI2B1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P09486


Technical Details

Immunogen	Full length human recombinant protein of human SPARC (NP_003109) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:4000 IHC 1:150


Anti-SPARC Mouse Monoclonal Antibody [Clone ID: OTI2B1] (M00862) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SPARC (TA507361)



Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-SPARC mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-SPARC mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPARC Mouse Monoclonal Antibody [Clone ID: OTI2B1]

For Research Use Only. Not for use in diagnostic procedures.