

Anti-Glucokinase GCK Monoclonal Antibody

Catalog Number: M00884


Overview

Product Name	Anti-Glucokinase GCK Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Glucokinase GCK Monoclonal Antibody catalog # M00884. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 4G6
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P35557


Technical Details

Immunogen	Purified recombinant fragment of human GCK expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Glucokinase GCK Monoclonal Antibody (M00884) Images



Western blotting validation for Anti-Glucokinase GCK Monoclonal Antibody M00884 Western Blot (WB) analysis using GCK Monoclonal Antibody against HEK293 (1) and GCK-hlgGfC transfected HEK293 (2) cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



ELISA validation of GCK using Anti-Glucokinase GCK Monoclonal Antibody (M00884). ELISA analysis of GCK antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glucokinase GCK Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.