

## Anti-RENT1 / hUPF1 Rabbit Monoclonal Antibody

Catalog Number: M00900-2


### Overview

Product Name	Anti-RENT1 / hUPF1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RENT1 / hUPF1 Rabbit Monoclonal Antibody catalog # M00900-2. Tested in WB, IHC, ICC/IF, IP, FC applications. This antibody reacts with Human, Mouse.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Monoclonal 19U80
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92900

### Technical Details

Immunogen	A synthesized peptide derived from human RENT1 / hUPF1
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:            Boster Bio's internal QC testing used WB 1:1000-1:5000            IHC 1:50-1:200            ICC/IF 1:50-1:200            IP 1:50            FC 1:50</p> <p>For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a></p>

## Anti-RENT1 / hUPF1 Rabbit Monoclonal Antibody (M00900-2) Images



Western blot analysis of RENT1 / hUPF1 expression in SH-SY5Y cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RENT1 / hUPF1 Rabbit Monoclonal Antibody