

Anti-GDF8 / Myostatin Reference Antibody (landogrozumab)

Catalog Number: M00926-4


Overview

Product Name	Anti-GDF8 / Myostatin Reference Antibody (landogrozumab)
Reactive Species	Bovine, Chicken, Human, Mouse, Rat, Turkey
Description	Boster Bio Anti-GDF8 / Myostatin Reference Antibody (landogrozumab) (Catalog # M00926-4). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Bovine, Turkey, Human, Mouse, Chicken. Endotoxin: < 0.494EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O14793


Technical Details

Isotype	IgG4
Form	Liquid

Anti-GDF8 / Myostatin Reference Antibody (landogrozumab) (M00926-4) Images



Anti-GDF8/Myostatin Reference Antibody (landogrozumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-GDF8/Myostatin Reference Antibody (landogrozumab) is more than 100%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GDF8 / Myostatin Reference Antibody (landogrozumab)

For Research Use Only. Not for use in diagnostic procedures.