

Anti-Tenascin C TNC Rabbit Monoclonal Antibody

Catalog Number: M00936-1

About TNC

Catalyzes the first step in leukotriene biosynthesis, and thereby plays a role in inflammatory processes.


Overview

Product Name	Anti-Tenascin C TNC Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Tenascin C TNC Rabbit Monoclonal Antibody catalog # M00936-1. Tested in WB application. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal CEG-20
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P24821

Technical Details

Immunogen	A synthesized peptide derived from human Tenascin C
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB 1:500-1:2000</p>

Anti-Tenascin C TNC Rabbit Monoclonal Antibody (M00936-1) Images



Western blot analysis of Cytotactin C expression in
(1)Human fetal brain lysate;(2)Human fetal kidney lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tenascin C TNC Rabbit Monoclonal Antibody