

Anti-PAX8 Antibody (Center)

Catalog Number: M00943-1

About PAX8

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.


Overview

Product Name	Anti-PAX8 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-PAX8 Antibody (Center) (Catalog # M00943-1). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 2131CT616.27.1
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q06710

Technical Details

Immunogen	This PAX8 antibody is generated from a mouse immunized with a recombinant protein from human PAX8.
Predicted Reactive Species	Human
Isotype	IgG1,kappa
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:4000

Anti-PAX8 Antibody (Center) (M00943-1) Images



All lanes : Anti-PAX8 Antibody (Center) at 1:4000 dilution
Lane 1: HeLa whole cell lysate
Lane 2: Raji whole cell lysate
Lane 3: SK-OV-3 whole cell lysate
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 48 kDa
Blocking/Dilution buffer: 5% NFD/MTBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAX8 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.