

## Anti-PAX-8 Rabbit Monoclonal Antibody

Catalog Number: M00943-4

### Overview

Product Name	Anti-PAX-8 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of PAX-8 protein. Validated for research with IHC,ELISA applications. Anti PAX-8 is reactive for Human samples
Application	ELISA, IHC
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	Q06710

### Technical Details

Isotype	IgG, Kappa
Purification	Protein A

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAX-8 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.