

Anti-Hsp70 Rabbit Monoclonal Antibody

Catalog Number: M00949-3


Overview

Product Name	Anti-Hsp70 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Hsp70 protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Hsp70 is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P0DMV8;P0DMV9


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:400-1:1000 WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-Hsp70 Rabbit Monoclonal Antibody (M00949-3) Images



Mouse brain was stained with Anti-Hsp70 rabbit antibody



Human prostate was stained with Anti-Hsp70 rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Hsp70 antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: A172 Lane 2: HeLa Lane 3: Mouse brain Lane 4: Rat brain
Predicted band size: 70kDa Observed band size: 70kDa

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hsp70 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.