

Anti-IRF5 Antibody (Monoclonal, 31118)

Catalog Number: M00958-4

About IRF5

Interferon regulatory factor 5, also called IRF5 or SLEB10, is a protein that in humans is encoded by the IRF5 gene. IRF5 gene is mapped to 7q32.1. This gene encodes a member of the interferon regulatory factor (IRF) family, a group of transcription factors with diverse roles, including virus-mediated activation of interferon, and modulation of cell growth, differentiation, apoptosis, and immune system activity. Members of the IRF family are characterized by a conserved N-terminal DNA-binding domain containing tryptophan (W) repeats. Multiple transcript variants encoding different isoforms have been found for this gene, and a 30-nt indel polymorphism (SNP rs60344245) can result in loss of a 10-aa segment. This gene is a transcription factor involved in the induction of interferons IFNA and INFB and inflammatory cytokines upon virus infection.

Overview

Product Name	Anti-IRF5 Antibody (Monoclonal, 31118)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IRF5 Antibody (Monoclonal, 31118) catalog # M00958-4. Tested in WB, IHC, IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, IHC, WB
Clonality	Monoclonal 31118
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13568

Technical Details

Immunogen	Recombinant protein within human IRF5 aa 3-198.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-200
Immunofluorescence, 1:50-200
ImmunoPrecipitation, 1:50

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IRF5 Antibody (Monoclonal, 31118)

For Research Use Only. Not for use in diagnostic procedures.