

## Anti-IL-22 Reference Antibody (fezakinumab)

Catalog Number: M00963-2


### Overview

Product Name	Anti-IL-22 Reference Antibody (fezakinumab)
Reactive Species	Cynomolgus monkey, Human, Rat
Description	Boster Bio Anti-IL-22 Reference Antibody (fezakinumab) (Catalog # M00963-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Cynomolgus monkey, Human. Endotoxin: < 0.499EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q9GZX6


### Technical Details

Isotype	IgG1
Form	Liquid


## Anti-IL-22 Reference Antibody (fezakinumab) (M00963-2) Images



Anti-IL-22 Reference Antibody (fezakinumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-22 Reference Antibody (fezakinumab) is more than 99.05%



Immobilized human IL 22

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-22 Reference Antibody (fezakinumab)

For Research Use Only. Not for use in diagnostic procedures.