

Anti-HDAC4/Histone Deacetylase 4 Rabbit Monoclonal Antibody

Catalog Number: M00971

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-HDAC4/Histone Deacetylase 4 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-HDAC4/Histone Deacetylase 4 Rabbit Monoclonal Antibody catalog # M00971. Tested in WB, ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Monoclonal HGI-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P56524

Technical Details

Immunogen	A synthesized peptide derived from human HDAC4
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

	ICC/IF 1:50-200
--	-----------------

Anti-HDAC4/Histone Deacetylase 4 Rabbit Monoclonal Antibody (M00971) Images



Western blot analysis of HDAC4 using anti-HDAC4 antibody (M00971). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Jurkat whole cell lysates, Lane 2: human THP-1 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-HDAC4 antigen affinity purified monoclonal antibody (M00971) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for HDAC4 at approximately 140 kDa. The expected band size for HDAC4 is at 119 kDa.

1 Publications Citing This Product

1. PubMed ID: 26473953, Arsenic trioxide reduces global histone H4 acetylation at lysine 16 through direct binding to histone acetyltransferase hMOF in human cells

Visit bosterbio.com/anti-hdac4-rabbit-monoclonal-antibody-m00971-boster.html to see all 1 publications.

Anti-HDAC4/Histone Deacetylase 4 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.