

Anti-TWIST1 monoclonal Antibody

Catalog Number: M00980


Overview

| | |
|----------------------|---|
| Product Name | Anti-TWIST1 monoclonal Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-TWIST1 monoclonal Antibody catalog # M00980. Tested in WB,IHC,ICC,Flow Cytometry applications. This antibody reacts with Human,Mouse. |
| Application | Flow Cytometry, IHC, ICC, WB |
| Clonality | Monoclonal 6H12 |
| Formulation | Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Mouse |
| Uniprot ID | Q15672 |


Technical Details

| | |
|---------------------|--|
| Immunogen | Purified recombinant fragment of human TWIST1 (AA: 9-74) expressed in E. Coli. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer |
| Suggested Dilutions | WB: 1:500-1:1000 IHC: 1:50-1:200 |

Anti-TWIST1 monoclonal Antibody (M00980) Images



Western blot (WB) analysis of TWIST1 monoclonal antibody
Lane1:NIH/3T3 whole cell lysate Lane2:Jurkat whole cell lysate Lane3:HeLa whole cell lysate Lane4: A549 whole cell lysate Lane5:Raji whole cell lysate Lane6:OCM1 whole cell lysate



Immunofluorescence analysis of HeLa cells using TWIST1 monoclonal antibody (green).

1 Publications Citing This Product

1. PubMed ID: 33237234, Rolim LSA, Mafra RP, Santos HBP, Souza LB, Pinto LP. Role of Twist and Podoplanin in Partial Epithelial-Mesenchymal Transition in Oral Squamous Cell Carcinoma. Braz Dent J. 2020 Nov-Dec;31(6):623-633. doi: 10.1590/0103-6440202003542. PMID: 33237234.

Visit bosterbio.com/anti-twist1-monoclonal-antibody-m00980-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TWIST1 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.