

Anti-PGP9.5 Rabbit Monoclonal Antibody

Catalog Number: M01018-1


Overview

Product Name	Anti-PGP9.5 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PGP9.5 Rabbit Monoclonal Antibody catalog # M01018-1. Tested in WB, IHC, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IHC, WB
Clonality	Monoclonal AOCE-21
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P09936


Technical Details

Immunogen	A synthesized peptide derived from human PGP9.5
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200 IP 1:50

Anti-PGP9.5 Rabbit Monoclonal Antibody (M01018-1) Images



Western blot analysis of PGP9.5 expression in SH-SY5Y cell lysate.



Immunohistochemical analysis of paraffin-embedded human brain, using PGP9.5 Antibody.

2 Publications Citing This Product

1. PubMed ID: 31305258, Age-associated changes of the intrinsic nervous system in relation with interstitial cells in the pre-weaning goat rumen Yu Liang,1 Imran Tarique,1 Waseem Ail Vistro,1 Yifei Liu,1 Ziyu Wang,2 Abdul haseeb,1 Noor Samad Gandahi,1 Adeela Iqbal,1 Siyi Wang,
2. PubMed ID: 15948235, Nude mice multi-drug resistance model of orthotopic transplantation of liver neoplasm and Tc-99m MIBI SPECT on p-glycoprotein

Visit bosterbio.com/anti-pgp9-5-rabbit-monoclonal-antibody-m01018-1-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PGP9.5 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.