

Anti-DCX/Doublecortin Rabbit Monoclonal Antibody

Catalog Number: M01053


Overview

Product Name	Anti-DCX/Doublecortin Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DCX/Doublecortin Rabbit Monoclonal Antibody catalog # M01053. Tested in WB, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal AOAD-4
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43602

Technical Details

Immunogen	A synthesized peptide derived from human DCX
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 FC 1:50-1:200

Anti-DCX/Doublecortin Rabbit Monoclonal Antibody (M01053) Images



Western blot analysis of DCX expression in C6 cell lysate.

2 Publications Citing This Product

1. PubMed ID: 20973978, Zhang Xm, Du F, Yang D, Yu Cj, Huang Xn, Liu W, Fu J. *Bmc Neurosci.* 2010 Oct 25;11:138. Doi: 10.1186/1471-2202-11-138. Transplanted Bone Marrow Stem Cells Relocate To Infarct Penumbra And Co-Express Endogenous Proliferative And Immature Neuronal M...
2. PubMed ID: 21699735, Zhang Xm, Du F, Yang D, Wang R, Yu Cj, Huang Xn, Hu Hy, Liu W, Fu J. *Bmc Neurosci.* 2011 Jun 24;12:61. Doi: 10.1186/1471-2202-12-61. Granulocyte Colony-Stimulating Factor Increases The Therapeutic Efficacy Of Bone Marrow Mononuclear Cell Transplant...

Visit bosterbio.com/anti-dcx-rabbit-monoclonal-antibody-m01053-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DCX/Doublecortin Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.