

Anti-RHD / CD240d Reference Antibody (LFB Anti-RhD)

Catalog Number: M01082-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-RHD / CD240d Reference Antibody (LFB Anti-RhD)
Reactive Species	Human
Description	Boster Bio Anti-RHD / CD240d Reference Antibody (LFB Anti-RhD) (Catalog # M01082-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.515EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q02161


Technical Details

Isotype	IgG1
Form	Liquid

Anti-RHD / CD240d Reference Antibody (LFB Anti-RhD) (M01082-1) Images



Anti-RHD/CD240d Reference Antibody (LFB Anti-RhD) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-RHD/CD240d Reference Antibody (LFB Anti-RhD) is more than 98.9%

Anti-RHD / CD240d Reference Antibody (LFB Anti-RhD)

For Research Use Only. Not for use in diagnostic procedures.