

Anti-Hsp90 alpha HSP90AA1 Rabbit Monoclonal Antibody

Catalog Number: M01103-3


Overview

Product Name	Anti-Hsp90 alpha HSP90AA1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Hsp90 alpha HSP90AA1 Rabbit Monoclonal Antibody catalog # M01103-3. Tested in WB, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, ICC, WB
Clonality	Monoclonal BAC-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07900


Technical Details

Immunogen	A synthesized peptide derived from human Hsp90 alpha
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200 IP 1:20

Anti-Hsp90 alpha HSP90AA1 Rabbit Monoclonal Antibody (M01103-3) Images



Western blot analysis of Hsp90 alpha expression in (1)PC-12 cell lysate; (2)HeLa cell lysate.



Immunofluorescent analysis of HeLa cells, using Hsp90 alpha Antibody .

2 Publications Citing This Product

1. PubMed ID: -, Youwen Zhang, Junye Wang, Shucheng Ye et al. Co-overexpression of STIP1 and Hsp90 correlates with progression and prognosis of lung adenocarcinoma, 29 January 2020, PREPRINT (Version 2) available at Research Square [https://doi.org/10.21203/rs.2.13680/v2]
2. PubMed ID: 21346199, Jiang Q, Wang Y, Li T, Shi K, Li Z, Ma Y, Li F, Luo H, Yang Y, Xu C. Mol Biol Cell. 2011 Apr 15;22(8):1167-80. Doi: 10.1091/Mbc.E10-10-0860. Epub 2011 Feb 23. Heat Shock Protein 90-Mediated Inactivation Of Nuclear Factor-??b Switches Autophagy To ...

Visit [bosterbio.com/anti-hsp90-alpha-rabbit-monoclonal-antibody-m01103-3-boster.html](https://www.bosterbio.com/anti-hsp90-alpha-rabbit-monoclonal-antibody-m01103-3-boster.html) to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hsp90 alpha HSP90AA1 Rabbit Monoclonal Antibody
For Research Use Only. Not for use in diagnostic procedures.