

Anti-RASSF1 Monoclonal Antibody

Catalog Number: M01104


Overview

Product Name	Anti-RASSF1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-RASSF1 Monoclonal Antibody catalog # M01104. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 6B9-D6-F8
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9NS23

Technical Details

Immunogen	Purified recombinant human RASSF1 protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

Anti-RASSF1 Monoclonal Antibody (M01104) Images



Western Blot analysis using RASSF1 Monoclonal Antibody against mouse liver, mouse brain lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RASSF1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.