

Anti-CD40L CD40LG Rabbit Monoclonal Antibody

Catalog Number: M01114

About CD40LG

Activated by icilin, eucalyptol, menthol, cold and modulation of intracellular pH. Involved in menthol sensation. Permeable for monovalent cations sodium, potassium, and cesium and divalent cation calcium.


Overview

Product Name	Anti-CD40L CD40LG Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CD40L CD40LG Rabbit Monoclonal Antibody catalog # M01114. Tested in WB application. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal AAIE-3
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P29965


Technical Details

Immunogen	A synthesized peptide derived from human CD40L
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB 1:500-1:2000


Anti-CD40L CD40LG Rabbit Monoclonal Antibody (M01114) Images




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of CD40L expression in MOLT-4 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD40L CD40LG Rabbit Monoclonal Antibody