

## Anti-TNFSF5 / CD40L / CD154 Reference Antibody (Toralizumab)

Catalog Number: M01114-7


### Overview

Product Name	Anti-TNFSF5 / CD40L / CD154 Reference Antibody (Toralizumab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TNFSF5 / CD40L / CD154 Reference Antibody (Toralizumab) (Catalog # M01114-7). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.533EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P29965


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-TNFSF5 / CD40L / CD154 Reference Antibody (Toralizumab) (M01114-7) Images



Anti-TNFSF5/CD40L/CD154 Reference Antibody (Toralizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFSF5/CD40L/CD154 Reference Antibody (Toralizumab) is more than 99.51%

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFSF5 / CD40L / CD154 Reference Antibody (Toralizumab)

For Research Use Only. Not for use in diagnostic procedures.