

## Anti-Podoplanin Rabbit Monoclonal Antibody

Catalog Number: M01124-4


### Overview

Product Name	Anti-Podoplanin Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Podoplanin protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Podoplanin is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	Q86YL7


### Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


## Anti-Podoplanin Rabbit Monoclonal Antibody (M01124-4) Images




Mosue lung was stained with anti-Podoplanin rabbit antibody



Rat lung was stained with anti-Podoplanin rabbit antibody



Human lung squamous carcinoma was stained with anti-Podoplanin rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Podoplanin antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: Mouse brain Lane 2: Rat brain Predicted band size: 17kDa Observed band size: 36kDa

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Podoplanin Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.