

## Anti-PBRM1 Antibody

Catalog Number: M01130-2

### About PBRM1

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Acts as a negative regulator of cell proliferation.


### Overview

Product Name	Anti-PBRM1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-PBRM1 Antibody (Catalog # M01130-2). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 1636CT387.122.52
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q86U86


### Technical Details

Immunogen	This PBRM1 antibody is generated from a mouse immunized with a recombinant protein of human PBRM1.
Predicted Reactive Species	Horse
Isotype	IgG1,k
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:1000

## Anti-PBRM1 Antibody (M01130-2) Images



Anti-PBRM1 Antibody at 1:500 dilution + HT-1080 whole cell lysate. Lysates/proteins at 20 µg per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 193 kDa. Blocking/Dilution buffer: 5% NFD/MTBST.



Anti-PBRM1 Antibody at 1:1000 dilution + HT-1080 whole cell lysate. Lysates/proteins at 20 µg per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 193 kDa. Blocking/Dilution buffer: 5% NFD/MTBST.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PBRM1 Antibody

For Research Use Only. Not for use in diagnostic procedures.