

Anti-Mps1 (TTK) Mouse Monoclonal Antibody [Clone ID: OTI1D4]

Catalog Number: M01132


Overview

Product Name	Anti-Mps1 (TTK) Mouse Monoclonal Antibody [Clone ID: OTI1D4]
Reactive Species	Human
Description	Boster Bio TTK mouse monoclonal antibody, clone OTI1D4. Catalog# M01132. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI1D4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P33981

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 189-464 of human TTK (NP_003309) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-Mps1 (TTK) Mouse Monoclonal Antibody [Clone ID: OTI1D4] (M01132) Images



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-TTK monoclonal antibody (1:500).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mps1 (TTK) Mouse Monoclonal Antibody [Clone ID: OTI1D4]

For Research Use Only. Not for use in diagnostic procedures.