

Anti-IL-alpha IL1A Monoclonal Antibody

Catalog Number: M01144


Overview

Product Name	Anti-IL-alpha IL1A Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL-alpha IL1A Monoclonal Antibody catalog # M01144. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1F3B3 / 3A6F11
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P01583

Technical Details

Immunogen	Purified recombinant fragment of human IL-1a expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-IL-alpha IL1A Monoclonal Antibody (M01144) Images



Western Blot analysis using IL-1alpha Monoclonal Antibody against truncated IL-1alpha recombinant protein.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-alpha IL1A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.