

## Anti-TNFRSF7 / CD27 Reference Antibody (Organon patent anti-CD27)

Catalog Number: M01148-5


### Overview

Product Name	Anti-TNFRSF7 / CD27 Reference Antibody (Organon patent anti-CD27)
Reactive Species	Human
Description	Boster Bio Anti-TNFRSF7 / CD27 Reference Antibody (Organon patent anti-CD27) (Catalog # M01148-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.539EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P26842


### Technical Details

Isotype	IgG2SA
Form	Liquid

## Anti-TNFRSF7 / CD27 Reference Antibody (Organon patent anti-CD27) (M01148-5) Images



Anti-TNFRSF7/CD27 Reference Antibody (Organon patent anti-CD27) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TNFRSF7/CD27 Reference Antibody (Organon patent anti-CD27) is more than 95%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF7 / CD27 Reference Antibody (Organon patent anti-CD27)

For Research Use Only. Not for use in diagnostic procedures.