

## Anti-IL-12b Reference Antibody (ebdarokimab)

Catalog Number: M01152-4

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




## Overview

|                      |                                                                                                                                                                                                                                    |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-IL-12b Reference Antibody (ebdarokimab)                                                                                                                                                                                       |
| Reactive Species     | Human                                                                                                                                                                                                                              |
| Description          | Boster Bio Anti-IL-12b Reference Antibody (ebdarokimab) (Catalog # M01152-4). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.540EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics                                                                                                                                                                                  |
| Clonality            | Monoclonal                                                                                                                                                                                                                         |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.                                                                                                    |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.                                                                                          |
| Uniprot ID           | P29460                                                                                                                                                                                                                             |


## Technical Details

|         |        |
|---------|--------|
| Isotype | IgG1   |
| Form    | Liquid |


## Anti-IL-12b Reference Antibody (ebdarokimab) (M01152-4) Images



Anti-IL-12b Reference Antibody (ebdarokimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-12b Reference Antibody (ebdarokimab) is more than 100%



Immobilized human IL-12b His at 2 μg/mL can bind Anti-IL-12b Reference Antibody (ebdarokimab)

### Anti-IL-12b Reference Antibody (ebdarokimab)

For Research Use Only. Not for use in diagnostic procedures.