

## Anti-GBA/Glucosylceramidase Rabbit Monoclonal Antibody

Catalog Number: M01162


### Overview

Product Name	Anti-GBA/Glucosylceramidase Rabbit Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-GBA/Glucosylceramidase Rabbit Monoclonal Antibody catalog # M01162. Tested in WB, IHC applications. This antibody reacts with Human, Rat.
Application	IHC, WB
Clonality	Monoclonal AACA-7
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04062


### Technical Details

Immunogen	A synthesized peptide derived from human GBA
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200


## Anti-GBA/Glucosylceramidase Rabbit Monoclonal Antibody (M01162) Images



All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.



Western blot analysis of GBA expression in U87-MG cell lysate.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GBA/Glucosylceramidase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.