

Anti-PKM Rabbit Monoclonal Antibody

Catalog Number: M01173-2


Overview

Product Name	Anti-PKM Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of PKM protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti PKM is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P14618


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-PKM Rabbit Monoclonal Antibody (M01173-2) Images




Rat lung was stained with Anti-PKM rabbit antibody




Mouse lung was stained with Anti-PKM rabbit antibody



Human kidney was stained with Anti-PKM rabbit antibody



Human lung was stained with Anti-PKM rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-PKM antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: A549 Lane 2: NIH3T3 Lane 3: C6 Predicted band size: 57kDa Observed band size: 57kDa

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKM Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.