

Anti-HGS Mouse Monoclonal Antibody [Clone ID: OTI3H7]

Catalog Number: M01174


Overview

Product Name	Anti-HGS Mouse Monoclonal Antibody [Clone ID: OTI3H7]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio HGS mouse monoclonal antibody, clone OTI3H7 (formerly 3H7). Catalog# M01174. Tested in IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, WB
Clonality	Monoclonal OTI3H7
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O14964


Technical Details

Immunogen	Full length human recombinant protein of human HGS (NP_004703) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:4000 IF 1:100

Anti-HGS Mouse Monoclonal Antibody [Clone ID: OTI3H7] (M01174) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HGS (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HGS.



Anti-HGS mouse monoclonal antibody (M01174) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HGS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HGS Mouse Monoclonal Antibody [Clone ID: OTI3H7]

For Research Use Only. Not for use in diagnostic procedures.