

Anti-SERCA2 ATP2A2 Rabbit Monoclonal Antibody

Catalog Number: M01176


Overview

Product Name	Anti-SERCA2 ATP2A2 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SERCA2 ATP2A2 Rabbit Monoclonal Antibody catalog # M01176. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal AOGE-1
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16615

Technical Details

Immunogen	A synthesized peptide derived from human SERCA2
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 FC 1:50

Anti-SERCA2 ATP2A2 Rabbit Monoclonal Antibody (M01176) Images



Western blot analysis of SERCA2 expression in HeLa cell lysate.

1 Publications Citing This Product

1. PubMed ID: 24324360, Guo Y, Zeng Qc, Zhang Cq, Zhang Xz, Li Rx, Wu Jm, Guan J, Liu L, Zhang Xc, Li Jy, Wan Zm. Int J Med Sci. 2013 Nov 5;10(13):1837-45. Doi: 10.7150/Ijms.6786. Ecollection 2013. Extracellular Matrix Of Mechanically Stretched Cardiac Fibroblasts Improv...

Visit bosterbio.com/anti-serca2-rabbit-monoclonal-antibody-m01176-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SERCA2 ATP2A2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.