

Anti-FOXL2 Rabbit Monoclonal Antibody

Catalog Number: M01185

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-FOXL2 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-FOXL2 Rabbit Monoclonal Antibody catalog # M01185. Tested in WB application. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 21F38
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P58012

Technical Details

Immunogen	A synthesized peptide derived from human FOXL2
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

Anti-FOXL2 Rabbit Monoclonal Antibody (M01185) Images



Western blot analysis of FOXL2 expression in Jurkat cell lysate.

1 Publications Citing This Product

1. PubMed ID: 10.1158/1078-0432.CCR-16-2180, HMGA2-FOXL2 Axis Regulates Metastases and Epithelial-to-Mesenchymal Transition of Chemoresistant Gastric Cancer

Visit bosterbio.com/anti-foxl2-rabbit-monoclonal-antibody-m01185-boster.html to see all 1 publications.

Anti-FOXL2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.