

Anti-Choline Acetyltransferase Rabbit Monoclonal Antibody

Catalog Number: M01192-2


Overview

Product Name	Anti-Choline Acetyltransferase Rabbit Monoclonal Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-Choline Acetyltransferase Rabbit Monoclonal Antibody catalog # M01192-2. Tested in WB, IHC applications. This antibody reacts with Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal 31C75
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P28329


Technical Details

Immunogen	A synthesized peptide derived from human Choline Acetyltransferase
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200

Anti-Choline Acetyltransferase Rabbit Monoclonal Antibody (M01192-2) Images



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Choline Acetyltransferase expression in Mouse brain cell lysate.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Choline Acetyltransferase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.