

Anti-TRX1 TXN Rabbit Monoclonal Antibody

Catalog Number: M01219

About TXN

Activated by icilin, eucalyptol, menthol, cold and modulation of intracellular pH. Involved in menthol sensation. Permeable for monovalent cations sodium, potassium, and cesium and divalent cation calcium.

Overview


Product Name	Anti-TRX1 TXN Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TRX1 TXN Rabbit Monoclonal Antibody catalog # M01219. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal AOGC-20
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10599

Technical Details

Immunogen	A synthesized peptide derived from human TRX1
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB 1:500-1:2000

IHC 1:50-1:200
ICC/IF 1:50-1:200
IP 1:50

Anti-TRX1 TXN Rabbit Monoclonal Antibody (M01219) Images



Western blot analysis of TRX1 expression in HepG2 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRX1 TXN Rabbit Monoclonal Antibody