

## Anti-ApoA-V APOA5 Monoclonal Antibody

Catalog Number: M01242


### Overview

Product Name	Anti-ApoA-V APOA5 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ApoA-V APOA5 Monoclonal Antibody catalog # M01242. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 2G1H11 / 1F1E8
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q6Q788


### Technical Details

Immunogen	Purified recombinant fragment of human ApoA-V expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-ApoA-V APOA5 Monoclonal Antibody (M01242) Images



Western blotting validation for Anti-ApoA-V APOA5 Monoclonal Antibody M01242 Western Blot (WB) analysis using ApoA-V Monoclonal Antibody against Human Serum (A) and ApoA-V recombinant protein (B). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



ELISA validation of APOA5 using Anti-ApoA-V APOA5 Monoclonal Antibody (M01242). ELISA analysis of ApoA-V antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ApoA-V APOA5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.