

Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1E3]

Catalog Number: M01242-1


Overview

Product Name	Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1E3]
Reactive Species	Human
Description	Boster Bio APOA5 mouse monoclonal antibody, clone OTI1E3 (formerly 1E3). Catalog# M01242-1. Tested in IF, WB. This antibody reacts with Human.
Application	IF, WB
Clonality	Monoclonal OTI1E3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q6Q788


Technical Details

Immunogen	Full length human recombinant protein of human APOA5 (NP_443200) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:4000


Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1E3] (M01242-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY APOA5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-APOA5.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-APOA5 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Immunofluorescent staining of HeLa cells using anti-APOA5 mouse monoclonal antibody (M01242-1).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1E3]

For Research Use Only. Not for use in diagnostic procedures.