

Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1A1]

Catalog Number: M01242-3


Overview

Product Name	Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1A1]
Reactive Species	Human
Description	Boster Bio APOA5 mouse monoclonal antibody, clone OTI1A1 (formerly 1A1). Catalog# M01242-3. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI1A1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q6Q788


Technical Details

Immunogen	Full length human recombinant protein of human APOA5 (NP_443200) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:4000 IHC 1:150


Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OTI1A1] (M01242-3) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY APOA5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-APOA5.



Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-APOA5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-APOA5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-APOA5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein A V (APOA5) Mouse Monoclonal Antibody [Clone ID: OT11A1]

For Research Use Only. Not for use in diagnostic procedures.