

## Anti-p57 Kip2 Rabbit Monoclonal Antibody

Catalog Number: M01244-2


### Overview

Product Name	Anti-p57 Kip2 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p57 Kip2 Rabbit Monoclonal Antibody catalog # M01244-2. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal 17C19
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49918


### Technical Details

Immunogen	A synthesized peptide derived from human p57 Kip2
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:20


## Anti-p57 Kip2 Rabbit Monoclonal Antibody (M01244-2) Images




Immunohistochemical analysis of paraffin-embedded Human thyroid cancer, using the Antibody at 1:50 dilution.



Western blot analysis on HeLa cell using p57 Kip2 Antibody.





Immunohistochemical analysis of paraffin-embedded Human ovarian cancer, using the Antibody at 1:50 dilution.



Immunohistochemical analysis of paraffin-embedded rat kidney, using p57 Kip2 Antibody .

Immunohistochemical analysis of paraffin-embedded Human tongue cancer, using the Antibody at 1:50 dilution.



Immunofluorescent analysis of HeLa cells treated with dexamethasone, using p57 Kip2 Antibody .

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p57 Kip2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.