

Anti-p57kip2 Rabbit Monoclonal Antibody

Catalog Number: M01244-3


Overview

Product Name	Anti-p57kip2 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of p57kip2 (ABT76R) rabbit mAb protein. Validated for research with WB,IHC,ELISA applications. Anti p57kip2 is reactive for Human samples
Application	ELISA, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P49918

Technical Details

Isotype	IgG, Kappa
Purification	Protein A

Anti-p57kip2 Rabbit Monoclonal Antibody (M01244-3) Images



Human kidney tissue was stained with anti-p57kip2 rabbit Antibody

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p57kip2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.