

Anti-Fetuin A Rabbit Monoclonal Antibody

Catalog Number: M01251-1


Overview

Product Name	Anti-Fetuin A Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fetuin A Rabbit Monoclonal Antibody catalog # M01251-1. Tested in WB, IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Monoclonal 26A17
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02765

Technical Details

Immunogen	A synthesized peptide derived from human Fetuin A
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:500-1:2000 IP 1:50</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-Fetuin A Rabbit Monoclonal Antibody (M01251-1) Images



Western blot analysis of Fetuin A expression in Human plasma lysate.

1 Publications Citing This Product

1. PubMed ID: 10.2147/OTT.S222426, Roles Of EAAT1, DHFR, And Fetuin-A In The Pathogenesis, Progression And Prognosis Of Chondrosarcoma

Visit bosterbio.com/anti-fetuin-a-rabbit-monoclonal-antibody-m01251-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fetuin A Rabbit Monoclonal Antibody