

Anti-Chromogranin A CHGA Rabbit Monoclonal Antibody

Catalog Number: M01256-1


Overview

Product Name	Anti-Chromogranin A CHGA Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Chromogranin A CHGA Rabbit Monoclonal Antibody catalog # M01256-1. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal AOIC-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10645

Technical Details

Immunogen	A synthesized peptide derived from human Chromogranin A
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 FC 1:20

Anti-Chromogranin A CHGA Rabbit Monoclonal Antibody (M01256-1) Images



Western blot analysis of Chromogranin A expression in SH-SY5Y cell lysate.

1 Publications Citing This Product

1. PubMed ID: 31732520, Glumac PM, Gallant JP, Shapovalova M, Li Y, Murugan P, Gupta S, Coleman IM, Nelson PS, Dehm SM, LeBeau AM. Exploitation of CD133 for the Targeted Imaging of Lethal Prostate Cancer. Clin Cancer Res. 2020 Mar 1;26(5):1054-1064. doi:10.1158/1078-0432.CCR-19-1659. Epub 20

Visit bosterbio.com/anti-chromogranin-a-rabbit-monoclonal-antibody-m01256-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Chromogranin A CHGA Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.