

## Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17)

Catalog Number: M01262-5


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17)  |
| Reactive Species     | Cynomolgus monkey, Human, Mouse, Rat  |
| Description          | Boster Bio Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17) (Catalog # M01262-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.557EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics   |
| Clonality            | Monoclonal  |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.   |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.   |
| Uniprot ID           | P40967  |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG1   |
| Form    | Liquid |

## Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17) (M01262-5) Images



Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17) is more than 96.2%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PMEL Reference Antibody (Novartis patent anti-PMEL17)

For Research Use Only. Not for use in diagnostic procedures.