

## Anti-Ku80 XRCC5 Rabbit Monoclonal Antibody

Catalog Number: M01275


### Overview

Product Name	Anti-Ku80 XRCC5 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Ku80 XRCC5 Rabbit Monoclonal Antibody catalog # M01275. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal FAE-24
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P13010


### Technical Details

Immunogen	A synthesized peptide derived from human Ku80
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:5000-1:10000 IHC 1:50-1:200 ICC/IF 1:100-1:500 IP 1:50


## Anti-Ku80 XRCC5 Rabbit Monoclonal Antibody (M01275) Images




All lanes use the Antibody at 1:2W dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded human colon, using Ku80 Antibody.



All lanes use the Antibody at 1:2W dilution for 1 hour at room temperature.



Western blot analysis of Ku80 expression in HeLa cell lysate.

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



**Anti-Ku80 XRCC5 Rabbit Monoclonal Antibody**

For Research Use Only. Not for use in diagnostic procedures.