

## Anti-Retinoid X Receptor alpha/RXRA Antibody Picoband® (monoclonal, 5E7) FITC Conjugated

Catalog Number: M01299-FITC

### About RXRA

Retinoid X receptor alpha (RXR-alpha), also known as NR2B1 (nuclear receptor subfamily 2, group B, member 1) is a nuclear receptor that in humans is encoded by the RXRA gene. Retinoid X receptors (RXRs) and retinoic acid receptors (RARs) are nuclear receptors that mediate the biological effects of retinoids by their involvement in retinoic acid-mediated gene activation. These receptors function as transcription factors by binding as homodimers or heterodimers to specific sequences in the promoters of target genes. The protein encoded by this gene is a member of the steroid and thyroid hormone receptor superfamily of transcriptional regulators. Alternative splicing of this gene results in multiple transcript variants.

### Overview

Product Name	Anti-Retinoid X Receptor alpha/RXRA Antibody Picoband® (monoclonal, 5E7) FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Monoclonal 5E7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P19793

### Technical Details

Immunogen	E. coli-derived human RXRA recombinant protein (Position: A226-T462).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG1
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 10.1539/joh.L8084, Effects of Chronic Noise Exposure on Spatial Learning and Memory of Rats in Relation to Neurotransmitters and NMDAR2B Alteration in the Hippocampus

Visit [bosterbio.com/anti-retinoid-x-receptor-alpha-rxra-picoband-trade-antibody-m01299-boster.html](http://bosterbio.com/anti-retinoid-x-receptor-alpha-rxra-picoband-trade-antibody-m01299-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Retinoid X Receptor alpha/RXRA Antibody (monoclonal, 5E7) - FITC

For Research Use Only. Not for use in diagnostic procedures.