

## Anti-TIM1 HAVCR1 Monoclonal Antibody

Catalog Number: M01306


### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-TIM1 HAVCR1 Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Description          | Boster Bio Anti-TIM1 HAVCR1 Monoclonal Antibody catalog # M01306. Tested in WB application. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                                                           |
| Application          | WB                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality            | Monoclonal AFBC-8                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q96D42                                                                                                                                                                                                                                                                                                                                                                                                                                             |

### Technical Details

|                     |                                                |
|---------------------|------------------------------------------------|
| Immunogen           | A synthesized peptide derived from human TIM1. |
| Isotype             | Rabbit IgG                                     |
| Form                | Liquid                                         |
| Concentration       | 0.5mg/ml                                       |
| Purification        | Affinity-chromatography                        |
| Suggested Dilutions | WB 1:500-1:2000                                |

## Anti-TIM1 HAVCR1 Monoclonal Antibody (M01306) Images



Western blot analysis of TIM 1 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TIM1 HAVCR1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.