

Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab)

Catalog Number: M01324-3


Overview

Product Name	Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab) (Catalog # M01324-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.563EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O14786


Technical Details

Isotype	IgG1
Form	Liquid


Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab) (M01324-3) Images




Anti-NRP1/VEGF165R/CD304 Reference Antibody (vesencumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-NRP1/VEGF165R/CD304 Reference Antibody (vesencumab) is more than 100%



Immobilized human NRP1 His at 2 µg/mL can bind Anti-NRP1/VEGF165R/CD304 Reference Antibody (vesencumab)



Human NRP1 HEK293 cells were stained with Anti-NRP1/VEGF165R/CD304 Reference Antibody (vesencumab) and negative control protein respectively

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab)

For Research Use Only. Not for use in diagnostic procedures.